

Product Name :
Specimen Cutting Machines

Product Code :
TN331



Description :

Specimen Cutting Machines

Technical Specification :

Specimen Cutting Machines
Can be use metallic, polymer, etc
Specification:
RPM- 1500 max
Wheel size: 33 x 2 x 32 mm
Sample size (max): 30 mm x30 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003