

Product Name :
Specimen Mounting Press

Product Code :
NLAB-STELTLAB190008



Description :

Specimen Mounting Press

Technical Specification :

This Machine is Automatic Mounting Press for metallographic specimen. It is suitable for all Thermosetting Materials, because of the cooling system. It can automatically finish with the heating temperature, the holding time and acting force. The operator can change 4 specification patterns according to the difficult specimens. It can make 2 specimens at one time.

SPECIFICATION:

Pattern specimen : ?25mm, ?30mm, ?40mm, ?50mm

Power: 220V 50Hz

Maximum consumption: 1800W

Pressure: 0~2Mpa

Temperature: 0~300°C

Holding time: 0~99minutes and 99 seconds

Out line dimension: 615 x 510x 500mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003