

**Product Name :**  
Specimen Polishing Machine

**Product Code :**  
NLAB-SCIENTIFICEQ28025



**Description :**

Specimen Polishing Machine

**Technical Specification :**

Specimen Polishing Machine Is Designed And Manufactured On The Basis Of International  
Advanced Metallographic Sample Preparation Methods And Sample Preparation Conditions

Two Working Disc. Facility

Simple Operation Makes It Very Easy To Handle

**Specifications**

Work Disc. : 2

Disc. Diameter : 200 Mm

Rotate Speed : 150 Rpm

Power Supply : 220 ~ 380V  $\pm 10\%$ , 50/60Hz

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003