

Product Name :
Specimen Polishing Machine

Product Code :
NLAB-SCIENTIFICEQ28025



Description :

Specimen Polishing Machine

Technical Specification :

Specimen Polishing Machine Is Designed And Manufactured On The Basis Of International
Advanced Metallographic Sample Preparation Methods And Sample Preparation Conditions

Two Working Disc. Facility

Simple Operation Makes It Very Easy To Handle

Specifications

Work Disc. : 2

Disc. Diameter : 200 Mm

Rotate Speed : 150 Rpm

Power Supply : 220 ~ 380V \pm 10%, 50/60Hz

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003