

Product Name :
Spectrometer Prism

Product Code :
NLAB-SCIENTIFICEQ25006



Description :

Spectrometer Prism

Technical Specification :

- Optically worked for spectrometer, equilateral prism with two faces polished.

Available Size:

- 25x25
- 32x32
- 38x38
- Crown Glass Prism: $n=1.51$
- Dense Flint Glass Prism: $n=1.62$
- Extra Flint Glass Prism: $n=1.6$

- Double Extra Dense Flint Glass Prism: $n=1.71-1.74$

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003