

Product Name :
Spectrometer (Simple)

Product Code :
NLAB-SCIENTIFICEQ25005



Description :

Spectrometer (Simple)

Technical Specification :

- Used primarily for the precise measurement of radius of curvature of curved surfaces such as lenses and mirrors.
- Stainless steel disc of about 40mm diameter, graduated 100 divisions, rests on the center on the center of a three pronged brass table, through a steel screw.
- Each prong carries a pointed leg of MS below it equidistant from each other.
- A vertical scale, graduated 10-0-10mm mounted on one of the prongs.
- Screw pitch 1mm, LC 0.01mm
- Screw pitch 0.05mm, LC 0.005mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003