

Product Name :
Spectrophotometer

Product Code :
NLAB-SCIENTIFICIND13005



Description :

Spectrophotometer

Technical Specification :

UV/Visible Spectrophotometer,

Wavelength: range 198 – 1000nm,

resolution 1nm,

accuracy ± 2 nm,

repeatability ± 0.5 nm,

spectral bandwidth 5nm;

transmittance: 0 to 199.9%;

absorbance -0.300 to 2.500 A,

resolution 0.1%T, 0.001A,

accuracy ± 1 %T, V0.01

Abs at 1.000 Absorbance, stay light <0.5% at 340nm and 220nm;

concentration: range -300 to 9999,

resolution selectable 1/0.1/0.01/0.001,

calibration –blank with a single standard or factor, factor 0.001 to 10000, standard 0.001 to 1000; light source – Zenon lamp,

outputs USB, Analogue, Rs232, internal printer,

Power 24V, units: multiple. Complete with user manual, dozen cuvettes, and protective cover

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003