


Product Name :
 Spectrum Analyser 1050MHz

Product Code :
 SPA-0004



Description :

Spectrum Analyser 1050MHz

Technical Specification :

0.1% Frequency Response
±20KHz (frequency of Center Frequency +200KHz)
±0.1% (frequency Mark Sweep width +100KHz)
Frequency Display Resolution

Frequency Sweep 100KHz/div with 1,2 and step system and 0Hz/div(0 frequency sweep) accuracy : ±10%
Frequency Rate 20KHz/div to 400KHz Frequency Bandwidth
Video Filter 4KHz
Amplitude Range 10p to 10dB
Frequency Accuracy ±2dB (±3dB/10dB for every step) ±2dB
Amplitude Accuracy ±0.2dB (±0.1dB intermediate frequency width)
Resolution 800x600 (Range) of CRT Screen
Resolution ±2dB (of Vertical Display) 250MHz signal
Amplitude Accuracy ±10dB (±10dB ste p) error: +/-1dB/10Db
Maximum Input Level 0dBm (input attenuator set at 0dB) +20dBm (input attenuator set at 40dB)
Overload Frequency Stability

Input Socket Type
Output Impedance
Power Supply ±10%;50Hz/60Hz

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003