

**Product Name :**  
Spectrum Analyzers Testing

**Product Code :**  
SPA-0001



**Description :**

Spectrum Analyzers Testing

**Technical Specification :**

Accuracy of Mark Frequency :  $\pm 0.1\%$  (frequency sweep width +100KHz)

Frequency Display Resolution : 100KHz (4 1/2 bit LED display)

Frequency Sweep Width : 0.1MHz/div -100KHz/div with 1,2 and step system and 0Hz/div (0 frequency sweep)  
accuracy :  $\pm 10\%$

Intermediate Frequency Bandwidth : 20KHz and 400KHz

Video Filter : (when on) 4KHz

Amplitude Range : -107~+13dB

Reference Level : -27dB~+13dB (10dB for every step) +2dB

Average Noise Level : -97dB (20KHz intermediate frequency width)

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003