

Product Name :
Speedy Moisture Testers

Product Code :
NLAB-SOILELAB140033



Description :

Speedy Moisture Testers

Technical Specification :

The unit consists of a pressure vessel with clamp for sealing cap, capacity is 20 gms. Rubber sealing gasket, pressure gauge calibrated in percentage moisture content 0-25% x 0.25% on the wet weight basis, an electronic balance for weighing sample, a scoop for measuring carbide reagent, a bottle of reagent, one cleaning brush and a set of 4 steel balls for thorough mixing . Complete in highly polished wooden carrying case with handle.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003