

Product Name :
Sphygmomanometer Aneroid Round Shaped Mobile

Product Code :
NLABMEDICAL-M3230002

Description :

Sphygmomanometer Aneroid Round Shaped Mobile India, Sphygmomanometer Aneroid Round Shaped Mobile manufacturer India, suppliers, Medical Sphygmomanometer Aneroid Round Shaped Mobile India, Hospital Sphygmomanometer Aneroid Round Shaped Mobile India

Technical Specification :

Sphygmomanometer Aneroid Round Shaped - Mobile, ABS Round dial face adjustable for easy reading, stand on aluminium tube.

With 5 plastic castors, metal basket, standard latex bulb, inflation valve with spring, standard end valve, PVC tubing.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003