

Product Name :
Sphygmomanometer Digital

Product Code :
NLABMEDICAL-M2410001



Description :

Sphygmomanometer Digital

Technical Specification :

Five languages voice prompt

Voice volume can be adjustable

30 sets of memories

Average value for the latest 3 readings

Automatically inflate to proper pressure when necessary

Display of systolic pressure, diastolic pressure and pulse rate

simultaneously

Automatic power off after 3 minutes of no operation

Low battery indication

One-touch operation

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003