

Product Name :
Sphygmomanometer

Product Code :
NLABMEDICAL-M3230001



Description :

Sphygmomanometer India, Sphygmomanometer manufacturer India, suppliers, Medical Sphygmomanometer India, Hospital Sphygmomanometer India

Technical Specification :

Calibrated nylon cuff, stop pin manometer, 2-tube adult size bladder.

Standard latex bulb, Special vinyl zipper case, Inflation Valve with spring,

Standard end valve.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003