

Product Name :
Naugra Lab Spirometer

Product Code :
NLAB-LABORATORYEM21011

NAUGRA[®]



Description :

Naugra Lab Spirometer

Technical Specification :

Spirometer

Capacity 6 Litre.

Complete with chain / chord compensates counter balance to float through pulley.

Tube calibrated to denote volume.

Inlet and outlet tubes provided.

Complete with corrugated tube and mouth piece

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003