

**Product Name :**  
Spittoon Mug

**Product Code :**  
NLABMEDICAL-M2950002



**Description :**

Spittoon Mug

**Technical Specification :**

Offered comes in quality stainless steel construction finish and finds use in areas like hospitals as well as medical care centers for collection of sputum.

Made available in stainless steel 304 grade finish, the superior designing of these mugs also makes them easy to clean and maintain.

Some of its features include durable construction finish, smooth mirror polished surface to provide for neat clean finish, durable handle and lid opening mechanism support among others.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003