


Product Name :
Spray Tester Machines

Product Code :
NLAB-TESTINGLABM87001



Description :

Spray Tester Machines

Technical Specification :

Spray Tester consists of a stand for holding the spray device the button portion have a sample holder at 45° whose center lie at the center line of the spray device. spray device made of glass (or plastic) laboratory funnel which is connected to the spray nozzle by a rubber tube.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003