

Product Name :
Spring Balance (Acrylic Tubular Type)

Product Code :
NLAB-SCIENTIFICEQ24003



Description :

Spring Balance (Acrylic Tubular Type)

Technical Specification :

Plastic balance made of clear acrylic tube. Each range is color coded for easy identification, and has dual scale graduated in both grams and Newton. The spring mechanism and working of the balance is clearly visible with provision for zero adjustment incorporated at the top. With suspension ring and load hook.

Ranges available are:

- 100g-1g
- 200g-2g
- 250g-2.5n

- 500g-5g
- 1000g-10g
- 2000g-20g
- 250g-2.5N
- 500g-5N
- 1kg-10N
- 2kg-20N
- 3kg-30N
- 5kg-50N

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003