

Product Name :
SS Grill Chicken Maker

Product Code :
NLAB-VOCATIONALTVET27082

Description :

SS Grill Chicken Maker

Technical Specification :

Frequency: 60 Hz

Material: Stainless Steel

Voltage: 220 V

Temperature Range: 100-250 Degree C

Weight: 4.6 kg

Features:

Longer life

Excellent finish

Superb design

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003