

Product Name :
Stainless Steel Ultrasonic Cleaner

Product Code :
NLABMEDICAL-M3430001



Description :

Stainless Steel Ultrasonic Cleaner

Technical Specification :

Ultrasonic Cleaner Tank Dim.:20cm x 15cm x 7cm liter

Cap.: 1.8 Ltrs U/S

Power: 50 Watts U/S

Freq.: 33 khz Input:180

Vac: 270 Vac

Ensures smooth working of cleaners within the voltage variation of 180-270 vac

Half the weight & size of the conventional machines

Competitive rates efficient service backup of India

Models can be made as per your specific requirements

Heater ensures high efficiency in cleaning

Cleaners are equipped with imported transducers for smooth and long life working

Industries Applicable:- Hospitals, Nursing Homes Educational Institutes

Electronic Industry, Jewelry Industry, Hotel Industry, Pharmaceutical Industry & Automobile Industry

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003