

**Product Name :**  
Stand Alone LCD Screen Digital Microscope

**Product Code :**  
NLAB-EDUSCHOOLEQ1730104

**Description :**

Stand Alone LCD Screen Digital Microscope

**Technical Specification :**

Focus Range: manual focus from 10mm to 50mm.

Image Sensor: 5M Pixels.

TV output: Available (to any monitor with TV in).

Light Source: 8 white LEDs with adjustable brightness.

Digital zoom: up to 4x.

Optical zoom: 5x to 50x, 125x.

LCD Display: 3.5".

Card Slot: Support MicroSD Card slot up to 32GB (MicroSD Card is not included in the package).

Magnification: 10x to 200x, 500x.

Bundled software with measurement functions.

Easily captures microscopic pictures and videos.

Adjustable stage makes looking at different sized objects easier.

Video Resolution: 640 x 480 pixels (AVI file, Frame Rate: 30FPS).

Photo Resolution: 5MP, 3MP, 1.3MP, 640 x 480 pixels (JPG file).

Battery: Rechargeable lithium battery (recharged by AC charger or USB) provides up to 3 hours standalone working time.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003