

Product Name : Stand Alone LCD Screen Digital Microscope	Product Code : NLAB-EDUSCHOOLEQ1730104
Description : Stand Alone LCD Screen Digital Microscope	
Technical Specification : Focus Range: manual focus from 10mm to 50mm. Image Sensor: 5M Pixels. TV output: Available (to any monitor with TV in). Light Source: 8 white LEDs with adjustable brightness. Digital zoom: up to 4x. Optical zoom: 5x to 50x, 125x. LCD Display: 3.5". Card Slot: Support MicroSD Card slot up to 32GB (MicroSD Card is not included in the package). Magnification: 10x to 200x, 500x. Bundled software with measurement functions. Easily captures microscopic pictures and videos. Adjustable stage makes looking at different sized objects easier. Video Resolution: 640 x 480 pixels (AVI file, Frame Rate: 30FPS). Photo Resolution: 5MP, 3MP, 1.3MP, 640 x 480 pixels (JPG file). Battery: Rechargeable lithium battery (recharged by AC charger or USB) provides up to 3 hours standalone working time.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003