

Product Name :
Standard Penetration Test Apparatus

Product Code :
NLAB-TESTINGLABM14003

Description :

Standard Penetration Test Apparatus

Technical Specification :

Features of Standard Penetration Test Apparatus :

Standard Penetration Tester.

Portable Penetrometer is used to determine the frequency of layering penetration of foundation. It can be used To determine the allowable bearing capacity of normal cohesive soil and plain fill based on the tested hammer bash No. (N10) Suitable for substantial field working, is more economical and convenient than large type ones

Technical Specifications :

Penetration force : 900 ~ 1100 N

Resolution : 1 N

Penetration Needle : 100 mm², 500mm², 200mm²

Control system : Fully digital

Display system : LCD / LED

Power Supply : 220 ~ 240 V AC, 50-60Hz

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com