

| | |
|--|--|
| Product Name : Standard Penetration Test Kit | Product Code : NLAB-TESTINGLABM41002 |
|--|--|

| |
|--|
| |
|--|

Description :

Standard Penetration Test Kit

Technical Specification :

Standard Penetration Test Kit

Drive Weight, 65 kg, fitted with chain.

Guide Pipe Assembly, fitted with a driving head on one side having standard “A” drill rod thread & a cap on the second side. Drive weight falls freely through a height of 75mm

Tripod stand,6 meter high connected by a tie bolt. Each leg is made of three part for easy transport. Complete with pulleys and stairing arrangements.

Tripod stand 5 meter long Each leg is made of two parts. Complete with pulley & stirring arrangement.

Pulley Spare.

“A” Drill Rod 1 meter long.

“A” Drill Rod 1.5 meter long.

“A” Drill Rod 3 meter long.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003