

**Product Name :**  
Standard Tar Viscometer

**Product Code :**  
NLAB-BTMENLAB150035



**Description :**

Standard Tar Viscometer

**Technical Specification :**

For determining the viscosity of cut back bitumen and road oil. The viscometer consists of a chrome plated copper bath, with a drain valve and a central tube to receive the test cup and to position the stirrer, and is mounted on a stand with levelling feet. Stirrer has a curved shield and is provided with an insulated handle, thermometer socket and swivel support for the valve.

Required Cup Sizes According to Different Standards:

EN 12846 : 2, 4 and 10 mm dia.

EN 13357 : 4 and 10 mm dia.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003