

**Product Name :**  
Stands Bossheads and Clamps

**Product Code :**  
NLAB-LBEMALAB230008



**Description :**

Stands Bossheads and Clamps

**Technical Specification :**

Laboratory clamps, stands and supports are used in combination for holding and securing other items of apparatus such as glassware, tubing, thermometers, probes, columns, funnels and overhead stirrers.

Laboratory clamps consist of a shaft and bosshead secured with thumbscrews onto a supporting rod or scaffold, and adjustable closing jaws for gripping objects safely and securely.

Most are metal clamps for the laboratory, with chemical and heat resistance, although the adjusting screws can be plastic, and the jaws may be lined with cork, plastic or rubber for enhanced grip

General purpose plastic clamps

Retort clamps for round bottom-flasks

Glass joint clamps for ground glassware

Three-prong extension clamps for irregular or outsized objects

Chain clamps, screen support clamps, and other specialised clamps for burettes, condensers or reaction vessels

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003