

**Product Name :**  
Steel Alloy Hook Connector

**Product Code :**  
NLAB-SAFETYLABP20062



**Description :**

Steel Alloy Hook Connector

**Technical Specification :**

Strength: 25 kN.

Opening: 56 mm.

Steel alloy hook.

Small and compact Drop Forged.

Weight: 349gm.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003