

Product Name :
Steel Cylinder For Argon

Product Code :
MU141



Description :

Steel Cylinder For Argon

Technical Specification :

Steel Cylinder For Argon

Capacity: 10 L

Low splatter and good weld appearance.

Usable in all positions.

Also suitable for pulsed arc and robotic use.

They use the same regulator fittings as BOC, Air products and Air liquid

No bottle rental

No environmental charges

No fees to come and collect a cylinder.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003