

**Product Name :**  
Steel Measure 400 ml

**Product Code :**  
NLAB-CEMETLAB130024



**Description :**

Steel Measure 400 ml

**Technical Specification :**

Made of steel, 88.1 mm height x 76.2 mm internal diameter, it is calibrated to contain  $400 \pm 1$  ml of water at 23°C.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003