

**Product Name :**

Stereo Zoom Articulated Arm Microscope

**Product Code :**

NLAB-EDUSCHOOLEQ173010

**Description :**

Stereo Zoom Articulated Arm Microscope

**Technical Specification :**

Eyepieces: Includes widefield 10x/20mm high eyepoint eyepieces. Interpupillary adjustment from 50mm - 75mm. Dual diopters on eyetubes. Eyepieces accept a 24mm diameter reticle. 15x and 20x eyepieces are available below.

Focusing: Focusing holder provides smooth focusing.

Total Magnification: 7x - 45x.

Head: Your choice of binocular or trinocular body viewing head rotates 360, locks in place by tightening set screw.

Working Distance: 100mm

Stand: Heavy duty articulated arm stand with tilting focusing holder and heavy duty base. Post on articulated arm stand is 15" tall.

Objective Lenses: Zoom range of 0.7x - 4.5x.

Illumination: Your choice of no illumination, LED ring light with 60 LED lights (includes ring light adapter), or dual pipe LED gooseneck illuminator.

Zoom Ratio: 6.4:1

Base Dimensions: 330mm (13") x 330mm (13") x 38mm (1.5")

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

