

**Product Name :**  
Stereo Zoom Articulated Arm Microscope

**Product Code :**  
NLAB-EDUSCHOOLEQ173010

**Description :**

Stereo Zoom Articulated Arm Microscope

**Technical Specification :**

**Eyepieces:** Includes widefield 10x/20mm high eyepoint eyepieces. Interpupillary adjustment from 50mm - 75mm. Dual diopters on eyetubes. Eyepieces accept a 24mm diameter reticle. 15x and 20x eyepieces are available below.

**Focusing:** Focusing holder provides smooth focusing.

**Total Magnification:** 7x - 45x.

**Head:** Your choice of binocular or trinocular body viewing head rotates 360, locks in place by tightening set screw.

**Working Distance:** 100mm

**Stand:** Heavy duty articulated arm stand with tilting focusing holder and heavy duty base. Post on articulated arm stand is 15" tall.

**Objective Lenses:** Zoom range of 0.7x - 4.5x.

**Illumination:** Your choice of no illumination, LED ring light with 60 LED lights (includes ring light adapter), or dual pipe LED gooseneck illuminator.

**Zoom Ratio:** 6.4:1

**Base Dimensions:** 330mm (13") x 330mm (13") x 38mm (1.5")

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

