

Product Name :
Stereo Zoom Boom Microscope

Product Code :
NLAB-EDUSCHOOLEQ173009

Description :

Stereo Zoom Boom Microscope

Technical Specification :

Eyepieces: Includes widefield 10x/20mm high eyepoint eyepieces. Interpupillary adjustment from 50mm - 75mm. Dual diopters on both eyetubes. Eyepieces accept a 24mm diameter reticle. 15x and 20x eyepieces are available below.

Head: Your choice of binocular or trinocular body viewing head rotates 360°, locks in place by tightening set screw.

Focusing: Includes focusing holder with coarse focus adjustment.

Total Magnification: 7x - 45x.

Stand: Heavy duty boom stand with square base. Horizontal post is 20.9" long. Vertical post is 15" tall and an 8.9" extension post is available. Lock ring on vertical post. Base dimensions are 11.2" x 32". Weight is 48 lbs.

Working Distance: 100mm

Illumination: Your choice of no illumination, LED ring light with 60 LED lights (includes ring light adapter), or dual pipe LED gooseneck illuminator.

Objective Lenses: Zoom range of 0.7x - 4.5x.

Zoom Ratio: 6.4:1

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003