

**Product Name :**  
Stereo Zoom Microscope

**Product Code :**  
NLAB-EDUSCHOOLEQ173059

**Description :**

Stereo Zoom Microscope

**Technical Specification :**

Series models are highly professional, extremely robust stereo microscopes which are easy to use with zoom function which can be operated from both sides Due to its nature, column supports, as well as the integrated dimmable LED/halogen with reflected and transmitted light. Furthermore it offers precise resolution, extra large visual field and is completely colour true. These models are ideally suited for use in the quality control department, assembly and repair workstations in the electronics and semiconductor industries as well as in biological laboratories and in research.

Stage: Large stage with handrest

Objectives: Long working distance

Zoom: Continuous

Optical system: Greenough zoom lenses

Illumination: LED/halogen/fluorescent transmitted and reflected

Viewing Head: 45° inclined 360° rotatable

Resolution: 225 Lp/mm @ 45X

Microscopy Techniques: Polarizing, Darkfield, oblique illumination

Stand: With two step focussing, universal, customized

Eyepiece: Wide field with dioptic adjustment

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## Naugralabequipments