

Product Name :
Stereoscopic Zoom Binocular and Trinocular Microscope

Product Code :
NLAB-EDUSCHOOLEQ173068

Description :

Stereoscopic Zoom Binocular and Trinocular Microscope

Technical Specification :

Nosepiece Quadruple nosepiece with accurate centring and positive click stops.

Stage Built in graduated mechanical stage of 135x120mm.

Observation Head 45-degrees inclined Trinocular head rotatable through 360-degrees.

Objectives Achromatic 4X - 10X- 40X SL and 100X SL Oil immersion Eyepieces.

Illumination 6V-20W halogen lamp, 220/110Volts & mirror attachment for working in day light.

Low coaxial drive controls for easy and smooth scan of specimen slide over range of 50x75mm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003