

Product Name :
Sterile Garment Storage Cabinet

Product Code :
NLAB-LABORATORYEM30001



Description :

Sterile Garment Storage Cabinet

Technical Specification :

These specially designed cabinets provide a wash of sterile filtered, clean air through HEPA Filters having an efficiency of 99.97% down to .3 micron particle size. This cabinet protects dust and other particulates on clean room garments or other materials.

For use in the change Room for garment storage & inside the clean room for storage of materials.

Self standing cabinet provides down flow clean air to clean room garments.

FEATURES Can be designed in Horizontal air flow system also.

All visible exposed surfaces white laminated, giving high impact resistance with smart clean appearance.

Large interior capacity for handling upto 20 garments at a time. Special provision for perforated shelves (optional) for component storage.

Low power consumption suitable for operation on single phase supply. Any size available.

APPLICATIONS

Pharmaceuticals Industries & Laboratories, Hospital Electronic Components and Semiconductor Industry, Photographic Laboratories, Optical Instruments Industries.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003