

Product Name :
Stirrup Bending Machine

Product Code :
NLAB-STELTLAB190012



Description :

Stirrup Bending Machine

Technical Specification :

With a revelation to our reputation and modify image we have established in this domain, we are offering a premium quality collection of Bar Bending Machine GF 20 to our renowned customers. Their sturdy designs, flawless finish and hassle free operations make these machines a preferred choice of our customers spread all over the globe. Moreover, we assure that all these products are in accordance with the needs and necessities of our respected customers.

Bending Rebar Diameter: Upto 16 mm

Plate Diameter: 750x600x10 mm

Bending Speed: 3 – 15 r/min

Motor Power: 1.5 kW, Voltage : 415 V

Motor Speed: 1440 r/min

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003