

**Product Name :**  
Stirrup

**Product Code :**  
NLAB-EDUSCHOOLEQ1740136



**Description :**

Stirrup

**Technical Specification :**

Stirrup made of brass wire for suspension of magnets in electrostatic experiments. 50 mm long & 45 mm high. Stirrup Supports a bar magnet or a friction rod, allowing it to turn in a horizontal plane. Use to determine polarity. Hooks are 50mm apart. Without magnet or cord

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003