

Product Name :
Stirrup

Product Code :
NLAB-EDUSCHOOLEQ1740136



Description :

Stirrup

Technical Specification :

Stirrup made of brass wire for suspension of magnets in electrostatic experiments. 50 mm long & 45 mm high. Stirrup Supports a bar magnet or a friction rod, allowing it to turn in a horizontal plane. Use to determine polarity. Hooks are 50mm apart. Without magnet or cord

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003