

**Product Name :**  
Stop Watch Analog

**Product Code :**  
NLAB-EDUSCHOOLEQ1740138



**Description :**

Stop Watch Analog

**Technical Specification :**

Stopwatch with 0-60 second main dial graduated to 1/5 second. The small dial reads up to 60 minutes and is graduated every 1 minute. Anti magnetic movement in a robust plated case with an overall dia. of 54mm and depth of 14mm approx.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003