

Product Name :
Stop Watch Analog

Product Code :
NLAB-EDUSCHOOLEQ1740138



Description :

Stop Watch Analog

Technical Specification :

Stopwatch with 0-60 second main dial graduated to 1/5 second. The small dial reads up to 60 minutes and is graduated every 1 minute. Anti magnetic movement in a robust plated case with an overall dia. of 54mm and depth of 14mm approx.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003