

**Product Name :**  
Storage and Instrument Cabinet Glass Door

**Product Code :**  
NLABMEDICAL-M2570010



**Description :**

Storage and Instrument Cabinet Glass Door

**Technical Specification :**

Specification :-

Instrument Cabinet is made of Mild Steel & Pre Treated Epoxy Powder Coated.

It has 4 Metal Shelves.

It has 2 Glass Doors.

It is fitted with Lock Key and Handle.

We are Manufacturer, Exporters and Suppliers of Storage & Instrument Cabinet Glass Door at tender and export competitive prices from India.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

