


Product Name :
Storage Container for Test Tubes

Product Code :
NLAB-EDUCATIONTVET15008



Description :

Storage Container for Test Tubes

Technical Specification :

Dimension 1500x750x750 mm

Grid, overflow and drain valve

Temperature regulation

Power 1500 W, power supply 220 V / 50 Hz

rack

Cadmium steel cover

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

