

**Product Name :**  
Strain Gauge AS Pressure Transducer Trainer

**Product Code :**  
NLAB-ELECTRONICSAB110034



**Description :**

Strain Gauge AS Pressure Transducer Trainer

**Technical Specification :**

**SCOPE OF LEARNING:**

Study of Pressure Measurement using Pressure Transducer

**TECHNICAL SPECIFICATIONS:**

Digital Meters:

Digital Voltmeter : 0 to 10V

Power Supplies:

DC Regulated Power Supply +5V, 250mA

Operated on Mains power 230V, 50Hz +10%

Components are mounted on the panels are:

Strain Gauge Transducer With Mechanical Arrangement

Pressure Gauge : 0 to 150 psi

Zero Adjust and Adder With Potentiometer Control

**SALIENT FEATURES:**

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Fuse for Short Circuit protection

Instruction manual.

Connections are brought out through 2mm Colored Sockets.

Patch Cords 2mm.

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x10"

## **Naugralabequipments**

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003