


**Product Name :**  
Strip Spin 12 Centrifuge with Rotor

**Product Code :**  
NLAB-EDUSCHOOLEQ45028



**Description :**

Strip Spin 12 Centrifuge with Rotor

**Technical Specification :**

Strip Spin 12 Mini Centrifuge from Naugralab is the only mini centrifuge designed for centrifugation of 12 position PCR strips.

Open the lid, and the rotor quickly decelerates for removal of samples.

Simply close the lid and the unit quickly ramps up to 5500 rpm.

With a fixed speed that produces 2,000 xg, this centrifuge is perfect for quick spin downs.

5,500 rpm / 2,000 xg	

<b>Capacity</b> 4 x 8 position PCR strips 48 x 0.2ml PCR tubes	
<b>Dimensions (W x D x H)</b> 14 x 20 x 11.2 cm	
<b>Weight</b> : / 1 kg	
<b>Electrical</b> Voltage 100-240V, 50-60 Hz, 300W	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003