

Product Name : Student Isometric Lever	Product Code : NLAB-LABORATORYEM77027
Description :	
Student Isometric Lever	
Technical Specification :	
The holder carrying a length of spring steel wire and a flat steel lever with holes for attaching thread wire. Very adjustable. Also available for Frontal Writing Points.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003