

Product Name :
Student Physiograph

Product Code :
NLAB-EDUSCHOOLEQ28002



Description :

Student Physiograph

Technical Specification :

**Salient
Features:-**

Unique
system of pen motors, Modular easily interchangeable in a minute or two.

Pen
length standard 120mm for ink writingsystem improving upon the international
standard (100mm).

Shorter
pen length would mean greater arcdistortions which is synonymous with waveform
distortion and arc distortion.

Pen
length below 100min is not internationally accepted.

Frequency Response:-With 120min ink penlength
70Hz-3db.

With
shorter ink pens 75 or 80min length this response will be appreciably higher.

Chart Drive:- With completely silent gear
train fitted with finest imported synchronous motor free from any backlash.

Noisy
or large humming gear train or motor can introduce Bio-feedback signals.

Z
fold chart paper for easy back reference and recording,

Ripple
free power supply-with high stabilization and overload protection.

Modular
Time & Event channel (optional).

Console:-

The
Horizontal section houses the 10 speed Chart Drive, Galvanometers (Modular)
with stabilized ripple free power supply & the vertical section housing the
amplifiers and with ink pens and inkwell assembly.

Ten
Chart Speeds provided are -50, 25, 10, 5, 2, 1, 0.5, 0.25, 0.1, 0.05 mm/sec.

Transducers:-

Pressure
Transducer:- Strain Gauge Using STATHAM MODEL P23 a range -30 to +300 mm/ Hg.
sufficient over pressure protection.

Applications:- Arterial and
Venous pressure.

Volume Transducer
(plethysmography):- For low volume changes (strain gauge type).

Heart,
Smooth Muscle, Isotonic muscle Isometric Force Transducer.

Respiration
Transducer.

Photo-electric
Pulse Transducer.

Isotonic
Transducer

Accessories:-

EKG
Limb electrodes with straps with 4 nos -1 set

Chest
electrode-1 no.

EEG
electrode- 8 nos.

Pens
ink writing - 2 nos.

Pens
cradles - 2 nos.

Ink
- 1/2 lit.

Chart
Paper ink writing -2 nos.

EEG
Jelly tube -1 tube.

EEG
paste -1 jar.

3
pin Junction Box -1 no.

Earthing
Cord -1 no.

Couplers:-

Strain
gauge coupler

AC
coupler

Pulse
Respiration coupler Isotonic coupler.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003