

Product Name : Study of Converters Analog to Digital and Digital To Analog	Product Code : NLAB-ELECTRONICSAB220034
Description : Study of Converters Analog to Digital and Digital To Analog	
Technical Specification : SCOPE OF LEARNING: Varification of 8-Bit Analog to Digital Conversion Varification of 8-Bit Digital to Analog Conversion TECHNICAL SPECIFICATIONS: Power Supplies: DC Power Supply +5V,150mA DC Power Supply +2.5V,150mA Dual DC Power Supply +15V,150mA DC Power Supply for Reference Voltage 0-5V,150mA Operated on Mains power 230V, 50Hz +10% Digital Meter: Voltmeter 2V/20V DC Technical specifications: Eight independent buffered logic level indicators for High / Low status Eight independent Output indicators for High / Low status ADC 0804 IC given on board with block diagram	

DAC MC1408 IC given on board with block diagram

SALIENT FEATURES:

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Instruction manual.

Connections are brought out through 2mm Colored Sockets.

Patch Cords 2mm.

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x8"

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003