

Product Name : Study of Parity Generator and Checker	Product Code : NLAB-ELECTRONICSAB220023
Description :	
Study of Parity Generator and Checker	
Technical Specification :	
SCOPE OF LEARNING:	
Parity Generator	
Parity Checker	
TECHNICAL SPECIFICATIONS:	
Power Supplies:	
DC Power Supply +5V,150mA	
Operated on Mains power 230V, 50Hz +10%	
Technical specifications:	
Five independent logic level inputs to select High / Low TTL levels, each with a LED to indicate high / low status and termination.	
Two independent buffered logic level indicators for High / Low status	
Circuit using 7486 IC.	
SALIENT FEATURES:	
Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.	
Instruction manual.	
Connections are brought out through 2mm Colored Sockets.	
Patch Cords 2mm.	

The trainer is housed in ABS Plastic cabinet/Wooden.

Size of the trainer set 12"x8"

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003