

Product Name : Study of Unijunction Transistor U. J. T.	Product Code : NLAB-ELECTRONICSAB420030
Description : Study of Unijunction Transistor (U. J. T.)	
Technical Specification : Experimental Training Board has been designed specifically to study different properties of Unijunction Transistor. Practical experience on this board carries great educative value for Science and Engineering Students. Object: 01. To plot V-I characteristics of a given U.J.T. (Unijunction Transistor). 02. To use the given U.J.T. as a Relaxation Oscillator. Features: The board consists of the following built-in parts : 01. Two 0-30V DC at 50mA, continuously variable regulated Power Supplies. 02. Digital Voltmeter DC 3½ Digit having range of 0-200V. 03. Digital Current meter DC 3½ Digit having range of 0-20mA 04. Uni-junction Transistor. 05. Adequate no. of other electronic components. 06. Mains ON/OFF switch, Fuse and Jewel light. 07. The unit is operative on 230V ±10% at 50Hz AC Mains.	

08. Adequate no. of patch cords stackable from rear both ends 4mm spring loaded plug length ½ meter.

Other Apparatus Required:

Cathode Ray Oscilloscope 20MHz

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003