

Product Name : Suction Set With Standard Handle	Product Code : NLABMEDICAL-M3010001
---	---



Description :
Suction Set With Standard Handle
Technical Specification :
Suction set used for suctioning body fluid during surgical procedure on thoracic cavity or abdominal cavity with the help of suction machine.
The product used for clear the secretion.
Made from transparent, medical grade PVC.
Handle design with round and crown tip with smooth edge.
Transparent tip allows visual monitoring of fluids.
Connecting tube design with universal connector.
Disposable, Sterile, single use.

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments