

Product Name :
Super Deluxe O.T. Table

Product Code :
NLABMEDICAL-M2760004

Description :

Super Deluxe O.T. Table

Technical Specification :

Versatile Super Deluxe O.T. Table is superior for general surgical procedures.

Accurate positions can be achieved with easy mechanical operation.

Five section stainless steel top & heavy base and column enhances the durability and hygiene.

No gap, any leakage & high maintenance ensures durability & stability of the table.

Key Structure:

The five section Stainless Steel table top is equipped with removable, antistatic mattresses.

All the positions are realized by a mechanical driver.

The table is furnished with hard, robust stainless steel shell which provides the sense of dirt free & sanitized operations.

The various positions such as Up-Down movements, Trendelenburg Reverse Trendelenburg.

Lateral Tilts, Flex- Reflex, Chair Position as realized by remote control.

Leg section & head section are interchangeable.

Technical Constraints:

Top Dimension: L1830 mm x W497 mm

Height Adjustment: 762 mm - 987mm

Trandelenburg / Reverse: 30° / 25°

Lateral Tilt: 20° / 20°

Kidney Elevator: 150 mm

Back Rest (up / Down): 80° / 25°

Leg Rest (up / Down): 15° / 90°

Head Rest (up / Down): 20° / 60°

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003